



Medical Wing Foundation

Trade and maintenance of medical devices



MESH LOCK

Mesh lock base provides one of the strongest bond strength available for one piece bracket.



MICRO -ETCHED BASE

All orthox brackets are micro Etched on the base. Instron machine delivers 9 mega pascal bond strength to orthox brackets.



Medical Wing Foundation


Trade and maintenance of medical devices



ROTH Description Brackets

Direct Bond .018 & 022

M/D	3.51	3.51	3.51	3.51	4.2	4.2	3.51	3.51	3.51	3.51
TORQUE	-7°	-7°	+2	+8	+12°	+12	+8	-2	+7°	-7
ANGLE	0°	0°	9°	9°	5°	5°	5°	9°	0°	0°
WITH HOOK .018X.028	10,01,17H	10,01,17H	10,01,15H					10,01,14H	10,01,16H	10,01,16H
WITHOUT HOOK .018X.028	10,01,17	10,01,17	10,01,15	10,01,12	10,01,10	10,01,11	10,01,13	10,01,14	10,01,16	10,01,16
WITH HOOK .022X.028	10,01,27H	10,01,27H	10,01,25H					10,01,24H	10,01,26H	10,01,26H
WITHOUT HOOK .022X.028	10,01,27	10,01,27	10,01,25	10,01,22	10,01,20	10,01,21	10,01,23	10,01,24	10,01,26	10,01,26
TOOTH	BICUSPIDS	BICUSPIDS	CUSPIDS	LATERAL	CENTRAL	CENTRAL	LATERAL	CUSPIDS	BICUSPIDS	BICUSPIDS



M/D	3.51	3.51	3.51	2.93	2.93	2.93	2.93	3.51	3.51	3.51
TORQUE	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
ANGLE	0°	0°	7°	0°	0°	0°	0°	7°	0°	0°
WITH HOOK .018X.028	10,01,36H	10,01,35H	10,01,34H					10,01,31H	10,01,32H	10,01,33H
WITHOUT HOOK .018X.028	10,01,36	10,01,35	10,01,34	10,01,30	10,01,30	10,01,30	10,01,30	10,01,31	10,01,32	10,01,33
WITH HOOK .022X.028	10,01,46H	10,01,45H	10,01,44H					10,01,41H	10,01,42H	10,01,43H
WITHOUT HOOK .022X.028	10,01,46	10,01,45	10,01,44	10,01,40	10,01,40	10,01,40	10,01,40	10,01,41	10,01,42	10,01,43
TOOTH	BICUSPIDS	BICUSPIDS	CUSPIDS	LATERAL	CENTRAL	CENTRAL	LATERAL	CUSPIDS	BICUSPIDS	BICUSPIDS

Brackets (5-5)

Qty - 20nos (1 case)

0.18 x 0.28	0.18 x 0.28	0.22 x 0.28	0.22 x 0.28
Hooks 3,4,5	Hooks 3	Hooks 3,4,5	Hooks 3
10,01,100H	10,01,100	10,01,200H	10,01,200

ORTHO



Medical Wing Foundation

Trade and maintenance of medical devices



BUCCAL TUBES



SINGLE

ORTHOX upper molar non-convertible tubes
Weldable .018 & .022

Single Tube

Tooth	Torque	Offset	In/Out	Length	0.18 x 0.25	0.22 x 0.28	Qty
67	-14	+8	1,07	3,5	20,11,10	20,12,10	1
67	-14	+8	1,07	3,5	20,11,11	20,12,11	1



ORTHOX Lower Molar Non-convertible Tubes
Weldable .018 & .022

Single Tube

Tooth	Torque	Offset	In/Out	Length	0.18 x 0.25	0.22 x 0.28	Qty
67	-14	+8	1,07	3,5	20,17,10	20,18,10	1
67	-14	+8	1,07	3,5	20,17,11	20,18,11	1



ORTHO



Medical Wing Foundation

Trade and maintenance of medical devices



BANDS

Lower Right

7 6



Lower Left

6 7



Upper Right

7 6



Upper Left

6 7



Size	Qty	Code 7 6	Code 6 7	Code 7 6	Code 6 7
29	1	30.05.10	30.06.10	30.07.10	30.08.10
29.5	1	30.05.11	30.06.11	30.07.11	30.08.11
30	1	30.05.12	30.06.12	30.07.12	30.08.12
30.5	1	30.05.13	30.06.13	30.07.13	30.08.13
31	1	30.05.14	30.06.14	30.07.14	30.08.14
31.5	1	30.05.15	30.06.15	30.07.15	30.08.15
32	1	30.05.16	30.06.16	30.07.16	30.08.16
32.5	1	30.05.17	30.06.17	30.07.17	30.08.17
33	1	30.05.18	30.06.18	30.07.18	30.08.18
33.5	1	30.05.19	30.06.19	30.07.19	30.08.19
34	1	30.05.20	30.06.20	30.07.20	30.08.20
34.5	1	30.05.21	30.06.21	30.07.21	30.08.21
35	1	30.05.22	30.06.22	30.07.22	30.08.22
35.5	1	30.05.23	30.06.23	30.07.23	30.08.23
36	1	30.05.24	30.06.24	30.07.24	30.08.24
36.5	1	30.05.25	30.06.25	30.07.25	30.08.25
37	1	30.05.26	30.06.26	30.07.26	30.08.26
37.5	1	30.05.27	30.06.27	30.07.27	30.08.27
38	1	30.05.28	30.06.28	30.07.28	30.08.28
38.5	1	30.05.29	30.06.29	30.07.29	30.08.29
39	1	30.05.30	30.06.30	30.07.30	30.08.30
39.5	1	30.05.31	30.06.31	30.07.31	30.08.31
40	1	30.05.32	30.06.32	30.07.32	30.08.32
40.5	1	30.05.33	30.06.33	30.07.33	30.08.33
41	1	30.05.34	30.06.34	30.07.34	30.08.34
41.5	1	30.05.35	30.06.35	30.07.35	30.08.35
42	1	30.05.36	30.06.36	30.07.36	30.08.36

ORTHO



Medical Wing Foundation
Trade and maintenance of medical devices



BANDS

Molar Master Kit 30.10.05

-

Pre-Welded Master Kit 30.10.06












Size	29,5	30	30,5	31	31,5	32	32,5	33	33,5	34	34,5	35	35,5	36	36,5	37	37,5	38	38,5	39	39,5	40	40,5	41	41,5	42	TOT
UR 6	2	2	2	2	4	4	4	4	4	5	5	5	8	8	10	10	10	10	8	4	4	2	2	2	2	2	125
UL 6	2	2	2	2	4	4	4	4	4	5	5	5	8	8	10	10	10	10	8	4	4	2	2	2	2	2	125
UR 6	2	2	2	2	4	4	4	4	4	5	5	5	8	8	10	10	10	10	8	4	4	2	2	2	2	2	125
UL 6	2	2	2	2	4	4	4	4	4	5	5	5	8	8	10	10	10	10	8	4	4	2	2	2	2	2	125
																										500	

WIRES

PREFORMED - NITINOL (NiTi) ARCH WIRES

Arch Forms

Preformed Super Elastic NiTi Arch Wires

Size	Qty	NATURAL		OVAL	
		Upper code	Lower code	Upper code	Lower code
 .012"	10	40.01.01	40.02.01	40.03.01	40.04.01
 .014"	10	40.01.02	40.02.02	40.03.02	40.04.02
 .016"	10	40.01.03	40.02.03	40.03.03	40.04.03
 .018"	10	40.01.04	40.02.04	40.03.04	40.04.04
 .020"	10	40.01.05	40.02.05	40.03.05	40.04.05
 .016" x .016"	10	40.01.06	40.02.06	40.03.06	40.04.06
 .016" x .022"	10	40.01.07	40.02.07	40.03.07	40.04.07
 .017" x .025"	10	40.01.08	40.02.08	40.03.08	40.04.08
 .018" x .025"	10	40.01.09	40.02.09	40.03.09	40.04.09
 .019" x .025"	10	40.01.10	40.02.10	40.03.10	40.04.10
 .021" x .025"	10	40.01.11	40.02.11	40.03.11	40.04.11

round

square

rectangular

WIRES



Medical Wing Foundation

Trade and maintenance of medical devices



Preformed Super Elastic Reverse NiTi Arch Wires

Shape	Size	Qty	NATURAL		OVAL	
			Upper code	Lower code	Upper code	Lower code
round	○ .012"	10	40,07,01	40,08,01	40,09,01	40,10,01
	○ .014"	10	40,07,02	40,08,02	40,09,02	40,10,02
	○ .016"	10	40,07,03	40,08,03	40,09,03	40,10,03
	○ .018"	10	40,07,04	40,08,04	40,09,04	40,10,04
square	□ .016" x .016"	10	40,07,05	40,08,05	40,09,05	40,10,05
	□ .016" x .022"	10	40,07,06	40,08,06	40,09,06	40,10,06
	□ .017" x .025"	10	40,07,07	40,08,07	40,09,07	40,10,07
	□ .019" x .025"	10	40,07,08	40,08,08	40,09,08	40,10,08
rectangular						

WIRES

Power Chain

Elastomeric Chain

Color	Continuous	Short	Long
Black	50,01,01	50,02,01	50,03,01
Blue	50,01,02	50,02,02	50,03,02
Clear	50,01,03	50,02,03	50,03,03
Fuchsia	50,01,04	50,02,04	50,03,04
Purple	50,01,05	50,02,05	50,03,05
Metallic Green	50,01,06	50,02,06	50,03,06
Red	50,01,07	50,02,07	50,03,07
Silver	50,01,08	50,02,08	50,03,08
Teal	50,01,09	50,02,09	50,03,09



Yellow	50,01,10	50,02,10	50,03,10
Pink	50,01,11	50,02,11	50,03,11
Orange	50,01,12	50,02,12	50,03,12

ELASTOMERICS

Ligature Ties

Latex-Free Ligature Ties Moulded

Color Code	Std. Pack	Color	Code
Black	50,04,01	Red	50,04,08
Blue	50,04,02	Silver	50,04,09
Yellow	50,04,03	Teal	50,04,10
Metallic Green	50,04,04	Orange	50,04,11
Light Green	50,04,05	White	50,04,12
Light Pink	50,04,06	Clear	50,04,13
Purple	50,04,07	Fuchsia	50,04,14
		Tooth	50,04,15



Latex-Free Ligature Ties Loose

Color Code	Std. Pack	Color Code	Std. Pack	Color Code	Std. Pack
Gray	50,04,16	Multi Colour	50,04,17	Clear	50,04,18

ORTHO

