

Medical Wing Foundation

Trade and maintenance of medical devices

apixiaTM
DIGITAL IMAGING

Scanner

Small Size, Big Impact

The Apixia PSP Scanner packs tremendous imaging power into a small package. Tiny enough to fit into the most crowded dental operatories for handy chairside use, yet robust enough to be placed in a central location for use by the entire office, the Apixia PSP Scanner is the ideal product for any size dental practice. Imaging with PSP plates offers all the flexibility and patient comfort of film, combined with the immediacy and environmental benefits of digital radiography. It's truly the best of both worlds.

Images in Just a Few Easy Steps

Imaging is fast and simple: Position and expose the flexible PSP plate just like a piece of film, remove it from the hygienic envelope, and drop it into the Apixia PSP Scanner. The smart top-loading design harnesses the power of gravity for simple loading and fuss-free use. In less than ten seconds, the crystal clear radiograph appears on the computer screen, ready for viewing. Then, just simply insert the next exposed plate when ready. The unit's in-line eraser clears the imaging plate for its next use? automatically! When the plate is ejected from the bottom of the machine, it's ready for immediate re-use. Plus, plates are reusable thousands of times to maximize your investment.



APIXIA PSP SCANNER

Dimensions	258mm x 138 mm x 183mm
Weight	3.6 kg
Voltage Requirements:	AC 100-240 V, AC 60/50HZ
Scan Time:	10 seconds
Resolution:	17 lp/mm
USB Cable Length:	1.8 meters
Connection Type:	USB 2.0
Computer System Requirements:	Windows XP, Windows Vista, Windows 7 32/64 bit
Note: <i>Computer not included</i>	512 MB, 1 GB recommended
Windows Versions	Pentium4 CPU or better
Min. RAM	SVGA 1024 x 768
Min. Video Card	40 GB
Min. Hard Drive Space	

X-Ray Scanner

