

Medical Wing Foundation

Trade and maintenance of medical devices

Digi-X
DENTAL IT UNBROKEN SYSTEM



■ SPECIFICATIONS

X-RAY Generator	Type Power Output	High Frequency Inverter 180W
X-RAY Installation	Tube Voltage Tube Current Heat Capacity	60kV 2mA 8.5Khu
X-RAY Tube	Anode Type Anode angle Focal Spot	Stationary 20° 0.8mm
Total Filtration		1.6mmAl
X-RAY Control	Exposure Time Set	0.02~0.18 sec (0.01 Step)
Power Requirement	Battery Battery Charger	DC 24V AC100~240V
Weight	Main Body Accessory	1.9Kg 0.6Kg
Dimension	Main body Charger	185 X 144 X 139 (mm) 100 X 65 X 35 (mm)

X-Ray

