

# Lotus Series

DISPOSABLE DENTAL HIGH-SPEED  
HANDPIECE WITH LED LIGHT

## LOTUS 401DB

产品描述：带灯一次性高速手机，3孔喷水，无需灭菌或维修，禁止循环使用

气压：180-250 KPa (1.8-2.5 kgf / cm<sup>2</sup>)

水压：198Kpa (2 kg)

转速： $\geq 300,000$  rpm

扭矩： $\geq 9.5$  N

夹持力：20-45N

插入插针力度： $\leq 45$ N

夹针方式：顶针式

适用车针： $\varnothing 1.59\text{mm}-\varnothing 1.6\text{mm} \times 21\text{mm}-23\text{mm}$  (直径×长度)

噪音： $\leq 70$ dB

光强度：20000 lx

光源：LED

灯泡电流：3.0-3.2 V

DESCRIPTION: With light disposable high speed handpiece Triple water spray  
Need no sterilization nor maintenance But circle use forbid

AIR PRESSURE: 180-250 KPa (1.8-2.5 kgf / cm<sup>2</sup>)

WATER PRESSURE: 198Kpa (2 kg)

ROTATION SPEED:  $\geq 300,000$  rpm

TORQUE:  $\geq 9.5$  N

GRIP POWER: 20-45N

POWER TO INSET BUR:  $\leq 45$ N

CHUCK TYPE: Friction Grip

BUR APPLICABLE:  $\varnothing 1.59\text{mm}-\varnothing 1.6\text{mm} \times 21\text{mm}-23\text{mm}$  (diameter×length)

NOISE:  $\leq 70$ dB

INTENSITY OF LIGHT: 20000 lx

SOURCE OF LIGHT: LED

CURRENT OF BULB: 3.0-3.2 V



## LOTUS 401D

产品描述：一次性高速手机，单孔喷水，无需灭菌或维修，禁止循环使用

气压：180-250 KPa (1.8-2.5 kgf / cm<sup>2</sup>)

水压：198Kpa (2 kg)

转速： $\geq 300,000$  rpm

扭矩： $\geq 9.5$  N

夹持力：20-45N

插入插针力度： $\leq 45$ N

夹针方式：顶针式

适用车针： $\varnothing 1.59\text{mm}-\varnothing 1.6\text{mm} \times 21\text{mm}-23\text{mm}$  (直径×长度)

噪音： $\leq 70$ dB

DESCRIPTION: Disposable high speed handpiece Single water spray Need

no sterilization nor maintenance But circle use forbid

AIR PRESSURE: 180-250 KPa (1.8-2.5 kgf / cm<sup>2</sup>)

WATER PRESSURE: 198Kpa (2 kg)

ROTATION SPEED:  $\geq 300,000$  rpm

TORQUE:  $\geq 9.5$  N

GRIP POWER: 20-45N

POWER TO INSET BUR:  $\leq 45$ N

CHUCK TYPE: Friction Grip

BUR APPLICABLE:  $\varnothing 1.59\text{mm}-\varnothing 1.6\text{mm} \times 21\text{mm}-23\text{mm}$  (diameter×length)

NOISE:  $\leq 70$ dB



BEING

Medical Equipment