

Medicine Wing Establishment

Trade and maintenance of medical devices



Low speed

ROSE202

DESCRIPTION	E type low speed handpiece 1:1 Ratio, with internal cooling system
AIR PRESSURE	245~392Kpa (2.5~4.0 kgf/cm ²)
WATER PRESSURE	198Kpa (2 kg)
ROTATION SPEED	22,000~27,000 rpm
CHUCK TYPE	Latch type
BUR APPLICABLE	Φ2.334~2.355 mm
NOISE	≤70dB



202CA(P)-B

202CA(P)

202SH-B

202SH

202AM

202AM-B

203EM

Dental Hand Pieces

