

*Medical Wing Foundation*

*Trade and maintenance of medical devices*

**YOUJOY®**

## Autoclave



physical space:	22L/17L
sterilizing temperature :	121°C, 134°C
Dry procedure :	vacuum drying
Display :	LCD,multi languages
test :	Bowie&Dick, Leak
sterilizing record :	There is a micro-printer optional, recording the whole process of sterilizing
Safety Features :	Pressure Release Valve and Computerized diagnosing and alarm system
water supply system :	Reservoirs(3.5L fresh water and 5 liter waste water reservoirs)
Voltage :	220V±10%
power :	2300W/1800W(10A)
instrument containers :	3 trays
Chamber size :	250mm*450mm/250mm*355mm
Net weight :	55Kg/52Kg

*Sterilization*



*Medical Wing Foundation*

*Trade and maintenance of medical devices*

**YOUJOY®**

## LUB90B

Hand piece lubricating device



Handpieces:3  
 Working pressure: 0.3-0.6Mpa  
 Voltage: 220V/50-60Hz  
 Air flow: 60L/min  
 Weight: 7.2kgs  
 Package size: 282x300x370mm

## SEAL-100

The best sealing scenario: Our SEAL 100 will change the sealing power to adjust it to voltage changes.  
 Stainless steel cover.  
 Two cutting devices.



Voltage	Power	Fuse TA	Sealing width	Press width	Weight	Packing size
V/Hz	W	TA	(max)mm	mm	kgs	mm
220/50-60	500	3	250	10	9.1	460*380*390

*Sterilization*



*Medical Wing Foundation**Trade and maintenance of medical devices***YOUJOY®****UC-01 Ultrasonic Cleaner**

Rinse time can be set from 2 to 30 min, rinse temperature can also be set. Power is automatically cut off when is lack of water.

External drainage system makes easy change of rinse liquid; low noise design makes quiet operation.

Specialized basket is designed for delicate minute devices that should be free from friction and braise, for example, high-speed and low-speed handpiece turbines.

Clear operation plate with digital display makes it easy to operate.

**Capacity:5L**

**Power:200W**

**Weight:7.7kgs**

**Packing size:410\*355\*300mm**

**DRINK**

Now with DRINK provided, you have a reliable source of high quality distilled water as long as you input pipe water isn't is simple.

Voltage: 220V

Frequency: 50/60HZ

Distiller water: 1.5L/H

Chamber size:  $\phi$  180x200

Outside size: 290x290x390 mm

Weight: 3.5kgs

Packing size: 245x285x510 mm

*Sterilization*